


FlexHD®

ACELLULAR HYDRATED DERMIS

Ready to perform breast reconstruction surgery?
Mentor is ready with FlexHD® Acellular Hydrated Dermis.



MENTOR—THE NAME YOU TRUST IN PLASTIC SURGERY—NOW BRINGS YOU A TRUSTED NAME IN ACELLULAR HUMAN DERMIS.

Through an alliance with MTF, the world's largest tissue bank, Mentor Worldwide LLC is proud to bring you **FlexHD®**. This acellular dermal matrix is derived from donated human allograft skin. Prehydrated and with no refrigeration necessary, **FlexHD®** is ready to use when you need it in your breast reconstruction surgeries.

CONSIDER THE BENEFITS OVER OTHER ACELLULAR HUMAN DERMIS:

- **Outstanding ease of use**
 - Hydrated and ready to use right out of the package
 - Reduces OR times and lowers operative costs
 - Provided sterile*
 - Exceptional handling and feel
- **Optimal results**
 - Demonstrates revascularization after 8 weeks^{1,2}
 - Available in matched pairs for improved bilateral reconstruction
 - When used in breast reconstruction, **FlexHD®** helps to define inframammary fold and provides an immediate breast pocket³

MTF Musculoskeletal
Transplant
Foundation
THE ALLOGRAFT LEADER™

*Passes USP <71> sterility tests.

PLEASE CONTACT YOUR MENTOR SALES
REPRESENTATIVE FOR MORE INFORMATION.

 **MENTOR®**
MAKE LIFE MORE *beautiful™*



ACELLULAR HYDRATED DERMIS

HOW TO PLACE AN ORDER FOR FlexHD®

Contact: MTF Customer Service Staff

- Telephone: 800-433-6576
- Fax: 732-661-2298
- Hours: Mon – Fri 8:15 AM – 7:30 PM (EST)

Tissue will be shipped next day by FedEx® at customer's expense.

HOW TO RETURN AN ORDER OF FlexHD®

Contact: MTF Customer Service Advocates

- Telephone: 800-946-9008 x 2307
- Hours: Mon – Fri 8:15 AM – 5:30 PM (EST)

All returns must be accompanied by a Return Authorization Number (RA#).

—RA#'s can be obtained by contacting MTF's Customer Service Advocates

—The twelve-digit serial number on the tissue packaging must be provided to obtain an RA#

Return requests must be made within 30 days of the original invoice date.

—Subject to a restocking charge of 20% of the service fee for the tissue

—Freight charges will not be credited

HOW TO REPORT A CUSTOMER COMPLAINT

Contact: MTF Customer Service Advocates

- Telephone: 800-946-9008 x 2307
- Hours: Mon – Fri 8:15 AM – 5:30 PM (EST)

1. Eberli D, Rodriguez S, Atala A, Yoo JJ. In vivo evaluation of acellular human dermis for abdominal wall repair. *J Biomed Mater Res A*. 2010;93(4):1527-1538. 2. Roth JS, Dexter DD, Lumpkins K, Bochicchio GV. Hydrated vs. freeze-dried human acellular dermal matrix for hernia repair: a comparison in a rabbit model. *Hernia*. 2009; 13(2):201-207. 3. Buck DW 2nd, Heyer K, DiBardino D, Bethke K, Kim JY. Acellular dermis-assisted breast reconstruction with the use of crescentric tissue expansion: a functional cosmetic analysis of 40 consecutive patients. *Aesthet Surg J*. 2010;30(2):194-200.

TISSUE CODE	DESCRIPTION	THICKNESS
ULTRA THIN		
479102	1 cm x 2 cm	0.2 mm - 0.4 mm
479204	2 cm x 4 cm	0.2 mm - 0.4 mm
479307	3 cm x 7 cm	0.2 mm - 0.4 mm
479407	4 cm x 7 cm	0.2 mm - 0.4 mm

THIN		
470102	1 cm x 2 cm	0.4 mm - 0.8 mm
470204	2 cm x 4 cm	0.4 mm - 0.8 mm
470307	3 cm x 7 cm	0.4 mm - 0.8 mm
470412	4 cm x 12 cm	0.4 mm - 0.8 mm

THICK		
471102	1 cm x 2 cm	0.8 mm - 1.7 mm
471104	1 cm x 4 cm	0.8 mm - 1.7 mm
471204	2 cm x 4 cm	0.8 mm - 1.7 mm
471207	2 cm x 7 cm	0.8 mm - 1.7 mm
471307	3 cm x 7 cm	0.8 mm - 1.7 mm
471510	5 cm x 10 cm	0.8 mm - 1.7 mm
471212	2 cm x 12 cm	0.8 mm - 1.7 mm
471312	3 cm x 12 cm	0.8 mm - 1.7 mm
471412	4 cm x 12 cm	0.8 mm - 1.7 mm
471512	5 cm x 12 cm	0.8 mm - 1.7 mm
471612	6 cm x 12 cm	0.8 mm - 1.7 mm
471416	4 cm x 16 cm	0.8 mm - 1.7 mm
471616	6 cm x 16 cm	0.8 mm - 1.7 mm
471816	8 cm x 16 cm	0.8 mm - 1.7 mm
471206	6 cm x 20 cm	0.8 mm - 1.7 mm
471208	8 cm x 20 cm	0.8 mm - 1.7 mm

ULTRA THICK		
472510	5 cm x 10 cm	≥1.8 mm
472312	3 cm x 12 cm	≥1.8 mm
472412	4 cm x 12 cm	≥1.8 mm
472512	5 cm x 12 cm	≥1.8 mm
472612	6 cm x 12 cm	≥1.8 mm
472416	4 cm x 16 cm	≥1.8 mm
472616	6 cm x 16 cm	≥1.8 mm
472816	8 cm x 16 cm	≥1.8 mm
472206	6 cm x 20 cm	≥1.8 mm
472208	8 cm x 20 cm	≥1.8 mm

BREAST KITS - 2 UNITS, THICK		
4B1612	6 cm x 12 cm	0.8 mm - 1.7 mm
4B1416	4 cm x 16 cm	0.8 mm - 1.7 mm
4B1616	6 cm x 16 cm	0.8 mm - 1.7 mm
4B1816	8 cm x 16 cm	0.8 mm - 1.7 mm

BREAST KITS - 2 UNITS, ULTRA THICK		
4B2612	6 cm x 12 cm	≥1.8 mm
4B2416	4 cm x 16 cm	≥1.8 mm
4B2616	6 cm x 16 cm	≥1.8 mm
4B2816	8 cm x 16 cm	≥1.8 mm